

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NSE2 Polyclonal Antibody

多克隆抗体

产品货号: PA05319

第三版

描述:FAM84B (Family With Sequence Similarity 84 Member B) is a Protein Coding gene. Diseases associated with FAM84B include Prostate Cancer Susceptibility. An important paralog of this gene is FAM84A.

抗原:Synthesized peptide derived from the N-terminal region of human NSE2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in esophageal squamous cell carcinomas.