

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat NSF Polyclonal Antibody

多克隆抗体

产品货号: PA05320

第三版

描述:NSF (N-Ethylmaleimide Sensitive Factor, Vesicle Fusing ATPase) is a Protein Coding gene. Diseases associated with NSF include Tetanus and Neuronal Intranuclear Inclusion Disease. Among its related pathways are Intra-Golgi traffic and Delta508-CFTR traffic / Sorting endosome formation in CF. GO annotations related to this gene include protein kinase binding and protein complex binding.

抗原:Synthesized peptide derived from the Internal region of human NSF

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: