

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NTPCR Polyclonal Antibody

多克隆抗体

产品货号: PA05323

第三版

描述:The protein encoded by this gene is a non-specific nucleoside triphosphatase that is slow-acting in vitro. This gene is overexpressed in many tumor tissues, and while it is not essential for the cell, overexpression is cytotoxic. However, the cytotoxicity is not related to its triphosphatase activity. NTPCR (Nucleoside-Triphosphatase, Cancer-Related) is a Protein Coding gene. Among its related pathways are Thiamine metabolism and superpathway of pyrimidine deoxyribonucleotides de novo biosynthesis. GO annotations related to this gene include poly(A) RNA binding and nucleoside-triphosphatase activity.

抗原:Synthesized peptide derived from the C-terminal region of human NTPase

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: