

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat/Monkey Nucleophosmin Polyclonal Antibody

多克隆抗体

产品货号: PA05325

第三版

**描述:** This gene encodes a phosphoprotein which moves between the nucleus and the cytoplasm. The gene product is thought to be involved in several processes including regulation of the ARF/p53 pathway. A number of genes are fusion partners have been characterized, in particular the anaplastic lymphoma kinase gene on chromosome 2. Mutations in this gene are associated with acute myeloid leukemia. More than a dozen pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants. NPM1 (Nucleophosmin) is a Protein Coding gene. Diseases associated with NPM1 include Leukemia, Acute Myeloid and Lymphomatoid Papulosis. Among its related pathways are Validated targets of C-MYC transcriptional activation and Chromosome Maintenance. GO annotations related to this gene include nucleic acid binding and protein homodimerization activity. An important paralog of this gene is NPM2.

**抗原:** Synthesized peptide derived from the N-terminal region of human Nucleophosmin

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) , 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:**