

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NUP160 Polyclonal Antibody

多克隆抗体

产品货号: PA05328

第三版

描述:NUP160 is 1 of up to 60 proteins that make up the 120-MD nuclear pore complex, which mediates nucleoplasmic transport. NUP160 (Nucleoporin 160) is a Protein Coding gene. Among its related pathways are HIV Life Cycle and Transport of the SLBP independent Mature mRNA. GO annotations related to this gene include nucleocytoplasmic transporter activity.

抗原:Synthesized peptide derived from the Internal region of human Nup160

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: