

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLCO1A2 Polyclonal Antibody

多克隆抗体

产品货号: PA05332

第三版

**描述:** This gene encodes a sodium-independent transporter which mediates cellular uptake of organic ions in the liver. Its substrates include bile acids, bromosulphophthalein, and some steroidal compounds. The protein is a member of the SLC21A family of solute carriers. Alternative splicing of this gene results in multiple transcript variants.

**抗原:** Synthesized peptide derived from the Internal region of human OATP1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**