

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Olfactory receptor 10X1 Polyclonal Antibody

多克隆抗体

产品货号: PA05362

第三版

**描述:**OR10X1 is an odorant receptor. Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. These proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes and share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals.

**抗原:**Synthesized peptide derived from the C-terminal region of human Olfactory receptor 10X1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4)，1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**