

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Olfactory receptor 13H1 Polyclonal Antibody

多克隆抗体

产品货号: PA05366

第三版

**描述:**OR13H1 belongs to the G protein coupled receptor 1 family and functions as an odorant receptor. Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G protein coupled receptors (GPCR) arising from single coding exon genes. Olfactory receptors share a 7 transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.

**抗原:**Synthesized peptide derived from the C-terminal region of human Olfactory receptor 13H1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**