

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Olfactory receptor 2AG1/2 Polyclonal Antibody

多克隆抗体

产品货号: PA05376

第三版

描述:OR2AG1 (Olfactory Receptor Family 2 Subfamily AG Member 1 (Gene/Pseudogene)) is a Protein Coding gene. Among its related pathways are Peptide ligand-binding receptors and Olfactory Signaling Pathway. GO annotations related to this gene include G-protein coupled receptor activity and olfactory receptor activity. An important paralog of this gene is OR2AG2. OR2AG2 (Olfactory Receptor Family 2 Subfamily AG Member 2) is a Protein Coding gene. Among its related pathways are Olfactory Signaling Pathway and Signaling by GPCR. GO annotations related to this gene include G-protein coupled receptor activity and olfactory receptor activity. An important paralog of this gene is OR2AG1.

抗原:Synthesized peptide derived from the Internal region of human Olfactory receptor 2AG1/2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: