

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Olfactory receptor 2D3 Polyclonal Antibody

多克隆抗体

产品货号: PA05379

第三版

描述:OR2D3 is an olfactory receptor. These receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G protein coupled receptors (GPCR) arising from single coding exon genes. Olfactory receptors share a 7 transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome.

抗原:Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2D3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: