

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat OSR2 Polyclonal Antibody

多克隆抗体

产品货号: PA05458

第三版

描述:OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors.OSR2 (Odd-Skipped Related Transcription Factor 2) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and sequence-specific DNA binding. An important paralog of this gene is OSR1.

抗原:Synthesized peptide derived from the Internal region of human OSR2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: