

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat p140 Cap Polyclonal Antibody

多克隆抗体

产品货号: PA05468

第三版

**描述:** SRCIN1 (SRC Kinase Signaling Inhibitor 1) is a Protein Coding gene. Diseases associated with SRCIN1 include Sclerosing Keratitis and Punctate Epithelial Keratoconjunctivitis. Among its related pathways are Adhesion and Neuroscience. GO annotations related to this gene include protein kinase binding and protein domain specific binding. An important paralog of this gene is KIAA1217.

**抗原:** Synthesized peptide derived from the Internal region of human p140Cap

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Expressed in some primary breast carcinomas where its presence is significantly associated with increased tumor size. Not detected in normal breast tissue.