

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human p27 Polyclonal Antibody

多克隆抗体

产品货号: PA05475

第三版

**描述:** Important regulator of cell cycle progression. Involved in G1 arrest. Potent inhibitor of cyclin E- and cyclin A-CDK2 complexes. Forms a complex with cyclin type D-CDK4 complexes and is involved in the assembly, stability, and modulation of CCND1-CDK4 complex activation. Acts either as an inhibitor or an activator of cyclin type D-CDK4 complexes depending on its phosphorylation state and/or stoichiometry.

**抗原:** Synthesized peptide derived from human p27 around the non-phosphorylation site of Thr198.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expressed in all tissues tested. Highest levels in skeletal muscle, lowest in liver and kidney.