

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PABPC3 Polyclonal Antibody

多克隆抗体

产品货号: PA05506

第三版

描述: Messenger RNA stability and translation initiation are extensively under the control of poly(A)-binding proteins (PABP). See PABPC1 (MIM 604679) for background information. PABPC3 (Poly(A) Binding Protein Cytoplasmic 3) is a Protein Coding gene. Among its related pathways are Deadenylation-dependent mRNA decay and RNA transport. GO annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is PABPC1.

抗原: Synthesized peptide derived from the Internal region of human PABP3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Testis specific.