

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SERPINA5 Polyclonal Antibody

多克隆抗体

产品货号: PA05509

第三版

描述:The protein encoded by this gene is a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. This gene is one in a cluster of serpin genes located on the q arm of chromosome 14. This family member is a glycoprotein that can inhibit several serine proteases, including protein C and various plasminogen activators and kallikreins, and it thus plays diverse roles in hemostasis and thrombosis in multiple organs.

抗原:Synthesized peptide derived from the Internal region of human PAI-3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed by the liver and secreted in plasma.