

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PARD3 Polyclonal Antibody

多克隆抗体

产品货号: PA05514

第三版

**描述:** This gene encodes a member of the PARD protein family. PARD family members interact with other PARD family members and other proteins; they affect asymmetrical cell division and direct polarized cell growth. Multiple alternatively spliced transcript variants have been described for this gene. PARD3 (Par-3 Family Cell Polarity Regulator) is a Protein Coding gene. Diseases associated with PARD3 include Gamma Heavy Chain Disease and Enamel Caries. Among its related pathways are Transcriptional activity of SMAD2/SMAD3-SMAD4 heterotrimer and Ras signaling pathway. GO annotations related to this gene include protein phosphatase binding and phosphatidylinositol-4,5-bisphosphate binding. An important paralog of this gene is PARD3B.

**抗原:** Synthesized peptide derived from the C-terminal region of human PARD3A

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) , 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Widely expressed.