

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse PAX2 Polyclonal Antibody

多克隆抗体

产品货号: PA05519

第三版

**描述:** PAX2 encodes paired box gene 2, one of many human homologues of the Drosophila melanogaster gene *prd*. The central feature of this transcription factor gene family is the conserved DNA-binding paired box domain. PAX2 is believed to be a target of transcriptional suppression by the tumor suppressor gene WT1. Mutations within PAX2 have been shown to result in optic nerve colobomas and renal hypoplasia. Alternative splicing of this gene results in multiple transcript variants.

**抗原:** Synthesized peptide derived from the Internal region of human Pax-2.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expressed in primitive cells of the kidney, ureter, eye, ear and central nervous system.