

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PLD4 Polyclonal Antibody

多克隆抗体

产品货号: PA05531

第三版

描述:PLD4 (Phospholipase D Family Member 4) is a Protein Coding gene. Among its related pathways are Role of phospholipids in phagocytosis and Innate Immune System. GO annotations related to this gene include phosphatidylinositol phospholipase C activity and N-acylphosphatidylethanolamine-specific phospholipase D activity. An important paralog of this gene is PLD3.

抗原:Synthesized peptide derived from the C-terminal region of human PC-PLD4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: