

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PDRG1 Polyclonal Antibody

多克隆抗体

产品货号: PA05543

第三版

描述:PDRG1 (P53 And DNA Damage Regulated 1) is a Protein Coding gene. GO annotations related to this gene include unfolded protein binding.

抗原:Synthesized peptide derived from the Internal region of human PDRG1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Predominantly expressed in normal testis and exhibits reduced but detectable expression in other organs.