

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PEX2 Polyclonal Antibody

多克隆抗体

产品货号: PA05553

第三版

**描述:** This gene encodes an integral peroxisomal membrane protein required for peroxisome biogenesis. The protein is thought to be involved in peroxisomal matrix protein import. Mutations in this gene result in one form of Zellweger syndrome and infantile Refsum disease. Alternative splicing results in multiple transcript variants encoding the same protein. PEX2 (Peroxisomal Biogenesis Factor 2) is a Protein Coding gene. Diseases associated with PEX2 include Peroxisome Biogenesis Disorder 5B and Peroxisome Biogenesis Disorder 5A. Among its related pathways are Peroxisome.

**抗原:** Synthesized peptide derived from the N-terminal region of human Peroxin 2

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**