

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PGLS Polyclonal Antibody

多克隆抗体

产品货号: PA05567

第三版

描述:PGLS (6-Phosphogluconolactonase) is a Protein Coding gene. Diseases associated with PGLS include Lymphatic System Disease. Among its related pathways are Pentose phosphate pathway and Pentose Phosphate Pathway (Erythrocyte). GO annotations related to this gene include monosaccharide binding and 6-phosphogluconolactonase activity.

抗原:Synthesized peptide derived from the C-terminal region of human PGLS

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: