

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ANP32B Polyclonal Antibody

多克隆抗体

产品货号: PA05569

第三版

**描述:**ANP32B (Acidic Nuclear Phosphoprotein 32 Family Member B) is a Protein Coding gene. GO annotations related to this gene include histone binding and RNA polymerase binding. An important paralog of this gene is ANP32A.

**抗原:**Synthesized peptide derived from the N-terminal region of human PHAPI2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in heart, lung, pancreas, prostate and in spleen, thymus and placenta.