
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse PHF3 Polyclonal Antibody

多克隆抗体

产品货号: PA05571

第三版

描述: This gene encodes a member of a PHD finger-containing gene family. This gene may function as a transcription factor and may be involved in glioblastomas development. Alternative splicing results in multiple transcript variants. PHF3 (PHD Finger Protein 3) is a Protein Coding gene. An important paralog of this gene is DIDO1.

抗原: Synthesized peptide derived from the C-terminal region of human PHF3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: