
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PHLDA3 Polyclonal Antibody

多克隆抗体

产品货号: PA05573

第三版

描述:PHLDA3 (Pleckstrin Homology Like Domain Family A Member 3) is a Protein Coding gene. Among its related pathways are Apoptosis and Autophagy and PI3K / Akt Signaling. GO annotations related to this gene include phosphatidylinositol-4,5-bisphosphate binding and phosphatidylinositol-3-phosphate binding. An important paralog of this gene is PHLDA1.

抗原:Synthesized peptide derived from the Internal region of human PHLDA3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed with lowest expression in liver and spleen.