

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse OLFML2B Polyclonal Antibody

多克隆抗体

产品货号: PA05575

第三版

**描述:**This gene encodes an olfactomedin domain-containing protein. Most olfactomedin domain-containing proteins are secreted glycoproteins. OLFML2B (Olfactomedin Like 2B) is a Protein Coding gene. GO annotations related to this gene include extracellular matrix binding. An important paralog of this gene is OLFML2A.

**抗原:**Synthesized peptide derived from the Internal region of human Olfactomedin-2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**