

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PIPOX Polyclonal Antibody

多克隆抗体

产品货号: PA05588

第三版

描述:PIPOX (Pipelicolic Acid And Sarcosine Oxidase) is a Protein Coding gene. Diseases associated with PIPOX include Peroxisomal Biogenesis Disorders. Among its related pathways are Lysine degradation and Peroxisome. GO annotations related to this gene include receptor binding and sarcosine oxidase activity.

抗原:Synthesized peptide derived from the Internal region of human PIPOX

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: