

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PITX1 Polyclonal Antibody

多克隆抗体

产品货号: PA05590

第三版

描述: This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. Members of this family are involved in organ development and left-right asymmetry. This protein acts as a transcriptional regulator involved in basal and hormone-regulated activity of prolactin. PITX1 (Paired Like Homeodomain 1) is a Protein Coding gene. Diseases associated with PITX1 include Clubfoot, Congenital, With Or Without Deficiency Of Long Bones And/Or Mirror-Image Polydactyly and Liebenberg Syndrome. Among its related pathways are DAG and IP3 signaling. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and transcriptional activator activity, RNA polymerase II core promoter proximal region sequence-specific binding. An important paralog of this gene is PITX2.

抗原: Synthesized peptide derived from the Internal region of human Pitx1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: