

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Monkey PXMP4 Polyclonal Antibody

多克隆抗体

产品货号: PA05605

第三版

描述:PXMP4 (Peroxisomal Membrane Protein 4) is a Protein Coding gene. Diseases associated with PXMP4 include Prostate Adenocarcinoma. Among its related pathways are Peroxisome.

抗原:Synthesized peptide derived from the N-terminal region of human Pmp24

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in normal prostate epithelial cells, and androgen-sensitive prostate adenocarcinoma cells. Not expressed in androgen-insensitive prostate adenocarcinoma cells.