

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Monkey POLR3A Polyclonal Antibody

多克隆抗体

产品货号: PA05609

第三版

描述:The protein encoded by this gene is the catalytic component of RNA polymerase III, which synthesizes small RNAs. The encoded protein also acts as a sensor to detect foreign DNA and trigger an innate immune response. POLR3A (RNA Polymerase III Subunit A) is a Protein Coding gene. Diseases associated with POLR3A include Leukodystrophy, Hypomyelinating, 7, With Or Without Oligodontia And/Or Hypogonadotropic Hypogonadism and Pol Iii-Related Leukodystrophies. Among its related pathways are Innate Immune System and RNA Polymerase III Transcription Initiation. GO annotations related to this gene include chromatin binding and RNA polymerase III activity.

抗原:Synthesized peptide derived from the N-terminal region of human POLR3A

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in the brain, in the cortex and the white matter (at protein level).