

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse POLR3E Polyclonal Antibody

多克隆抗体

产品货号: PA05611

---

第三版

**描述:** POLR3E (RNA Polymerase III Subunit E) is a Protein Coding gene. Among its related pathways are Innate Immune System and RNA Polymerase III Transcription Initiation. GO annotations related to this gene include DNA-directed 5-3 RNA polymerase activity and RNA polymerase III activity.

**抗原:** Synthesized peptide derived from the Internal region of human POLR3E

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**