

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat POLR3H Polyclonal Antibody

多克隆抗体

产品货号: PA05612

第三版

描述: POLR3H (RNA Polymerase III Subunit H) is a Protein Coding gene. Among its related pathways are Innate Immune System and RNA Polymerase III Transcription Initiation. GO annotations related to this gene include DNA-directed 5-3 RNA polymerase activity and RNA polymerase III activity.

抗原: Synthesized peptide derived from the C-terminal region of human POLR3H

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: