

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Monkey PPP4R2 Polyclonal Antibody

多克隆抗体

产品货号: PA05617

第三版

描述:The protein encoded by this gene is a regulatory subunit of the serine/threonine-protein phosphatase 4 complex. In addition to being required for efficient DNA double strand break repair, this complex plays a role in organization of microtubules at centrosomes and processing of spliceosomal snRNPs. Several transcript variants encoding different isoforms have been found for this gene. PPP4R2 (Protein Phosphatase 4 Regulatory Subunit 2) is a Protein Coding gene. Among its related pathways are DNA Double-Strand Break Repair. GO annotations related to this gene include protein binding, bridging.

抗原:Synthesized peptide derived from the Internal region of human PP4R2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed.