

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PPP1R14D Polyclonal Antibody

多克隆抗体

产品货号: PA05618

第三版

**描述:**Protein phosphatase-1 (PP1; see MIM 176875) is a major cellular phosphatase that reverses serine/threonine protein phosphorylation. PPP1R14D is a PP1 inhibitor that itself is regulated by phosphorylation. PPP1R14D (Protein Phosphatase 1 Regulatory Inhibitor Subunit 14D) is a Protein Coding gene. GO annotations related to this gene include protein phosphatase inhibitor activity. An important paralog of this gene is PPP1R14A.

**抗原:**Synthesized peptide derived from the Internal region of human PPP1R14D

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4)，1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Detected in colon, intestine, kidney and brain cortex.