

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse PKNOX2 Polyclonal Antibody

多克隆抗体

产品货号: PA05626

第三版

描述: Homeodomain proteins are sequence-specific transcription factors that share a highly conserved DNA-binding domain and play fundamental roles in cell proliferation, differentiation, and death. PKNOX2 belongs to the TALE (3-amino acid loop extension) class of homeodomain proteins characterized by a 3-amino acid extension between alpha helices 1 and 2 within the homeodomain. PKNOX2 (PBX/Knotted 1 Homeobox 2) is a Protein Coding gene. Diseases associated with PKNOX2 include Substance Dependence and Hydrolethalus Syndrome. GO annotations related to this gene include sequence-specific DNA binding and actin monomer binding. An important paralog of this gene is PKNOX1.

抗原: Synthesized peptide derived from the C-terminal region of human PREP-2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: