

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BTC Polyclonal Antibody

多克隆抗体

产品货号: PA05629

第三版

**描述:** This gene encodes a member of the epidermal growth factor (EGF) family of proteins. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the secreted growth factor. A secreted form and a membrane-anchored form of this protein bind to multiple different EGF receptors. This protein promotes pancreatic cell proliferation and insulin secretion, as well as retinal vascular permeability. Mutations in this gene may be associated with type 2 diabetes in human patients.

**抗原:** Synthesized peptide derived from the Internal region of human Probetacellulin

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Synthesized in several tissues and tumor cells. Predominantly expressed in pancreas and small intestine.