
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PRPF19 Polyclonal Antibody

多克隆抗体

产品货号: PA05631

第三版

描述: PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008 [PubMed 18263876]).

抗原: Synthesized peptide derived from the Internal region of human PRP19

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitous. Weakly expressed in senescent cells of different tissue origins. Highly expressed in tumor cell lines.