

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human PRSS33 Polyclonal Antibody

多克隆抗体

产品货号: PA05634

---

第三版

**描述:**PRSS33 (Protease, Serine 33) is a Protein Coding gene. GO annotations related to this gene include serine-type endopeptidase activity and serine-type peptidase activity. An important paralog of this gene is PRSS27.

**抗原:**Synthesized peptide derived from the Internal region of human PRSS33

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Predominantly expressed in macrophages. Present in the spleen, small and large intestine, lung and brain (at protein level). Highly expressed in peripheral leukocytes, ovary, retina, spleen and stomach. Moderately expressed in thymus, uterus and platelets, as well as some brain tissues, such as thalamus and fetal brain.