

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PRRX1 Polyclonal Antibody

多克隆抗体

产品货号: PA05636

第三版

描述:The DNA-associated protein encoded by this gene is a member of the paired family of homeobox proteins localized to the nucleus. The protein functions as a transcription co-activator, enhancing the DNA-binding activity of serum response factor, a protein required for the induction of genes by growth and differentiation factors. The protein regulates muscle creatine kinase, indicating a role in the establishment of diverse mesodermal muscle types. Alternative splicing yields two isoforms that differ in abundance and expression patterns.PRRX1 (Paired Related Homeobox 1) is a Protein Coding gene. Diseases associated with PRRX1 include Agnathia-Otocephaly Complex and Dysgnathia Complex. GO annotations related to this gene include sequence-specific DNA binding and HMG box domain binding. An important paralog of this gene is PRRX2.

抗原:Synthesized peptide derived from the Internal region of human PRX1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: