

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse PSKH1 Polyclonal Antibody

多克隆抗体

产品货号: PA05640

---

第三版

**描述:**PSKH1 (Protein Serine Kinase H1) is a Protein Coding gene. Among its related pathways are mRNA Splicing - Major Pathway. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is PSKH2.

**抗原:**Synthesized peptide derived from the Internal region of human PSK-H1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in all tissues and cell lines tested with the highest level of abundance in testis.