

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat RPL10 Polyclonal Antibody

多克隆抗体

产品货号: PA05661

第三版

描述: This gene encodes a ribosomal protein that is a component of the 60S ribosome subunit. The related protein in chicken can bind to c-Jun and can repress c-Jun-mediated transcriptional activation. Some studies have detected an association between variation in this gene and autism spectrum disorders, though others do not detect this relationship. There are multiple pseudogenes of this gene dispersed throughout the genome. Alternative splicing results in multiple transcript variants. RPL10 (Ribosomal Protein L10) is a Protein Coding gene. Diseases associated with RPL10 include Autism X-Linked 5 and X-Linked Intellectual Disability-Cerebellar Hypoplasia-Spondylo-Epiphyseal Dysplasia Syndrome. Among its related pathways are Viral mRNA Translation and rRNA processing in the nucleus and cytosol. GO annotations related to this gene include poly(A) RNA binding and structural constituent of ribosome. An important paralog of this gene is RPL10L.

抗原: Synthesized peptide derived from the N-terminal region of human QM.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: