

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RAB40B Polyclonal Antibody

多克隆抗体

产品货号: PA05671

第三版

描述:RAB40B is similar to a yeast protein which suggests a role of the gene product in regulating secretory vesicles. It may be a substrate recognition component of a SCF like ECS (Elongin Cullin SOCS box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

抗原:Synthesized peptide derived from the C-terminal region of human Rab 40B

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: