

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse RAB11FIP4 Polyclonal Antibody

多克隆抗体

产品货号: PA05677

---

第三版

**描述:**Proteins of the large Rab GTPase family (see RAB1A; MIM 179508) have regulatory roles in the formation, targeting, and fusion of intracellular transport vesicles. RAB11FIP4 is one of many proteins that interact with and regulate Rab GTPases (Hales et al., 2001 [PubMed 11495908]). The RAB11FIP4 effector protein acts in endosomal trafficking.

**抗原:**Synthesized peptide derived from the Internal region of human Rab11-FIP4

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Present at high level in testis (at protein level). Weakly expressed in other tissues.