
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat FRK Polyclonal Antibody

多克隆抗体

产品货号: PA05706

第三版

描述:Positively regulates PTEN protein stability through phosphorylation of PTEN on 'Tyr-336', which in turn prevents its ubiquitination and degradation possibly by reducing its binding to NEDD4. May function as a tumor suppressor.

抗原:Synthesized peptide derived from the Internal region of human Rak

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Restricted to cells lines derived from tissues of lymphoid, brain, breast, colon and bladder origin.