

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat RALY Polyclonal Antibody

多克隆抗体

产品货号: PA05707

第三版

描述:This gene encodes a member of the heterogeneous nuclear ribonucleoprotein (hnRNP) gene family. This protein may play a role in pre-mRNA splicing and in embryonic development. Alternate splicing results in multiple transcript variants. RALY (RALY Heterogeneous Nuclear Ribonucleoprotein) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and RNA binding. An important paralog of this gene is RALYL.

抗原:Synthesized peptide derived from the Internal region of human Raly

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in heart, brain, lung, liver, skeletal muscle, kidney and pancreas. Weakly expressed in placenta.