

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse RAP2C Polyclonal Antibody

多克隆抗体

产品货号: PA05711

第三版

描述:The protein encoded by this gene is a member of the Ras-related protein subfamily of the Ras GTPase superfamily. Members of this family are small GTPases that act as molecular switches to regulate cellular proliferation, differentiation, and apoptosis. This protein has been reported to activate in vitro transcriptional activity of the serum response element. Alternative splicing results in multiple transcript variants. RAP2C (RAP2C, Member Of RAS Oncogene Family) is a Protein Coding gene. Among its related pathways are Innate Immune System and ERK Signaling. GO annotations related to this gene include GTP binding and transcription coactivator activity. An important paralog of this gene is RAP2A.

抗原:Synthesized peptide derived from the Internal region of human Rap 2C

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in liver, skeletal muscle, prostate, uterus, rectum, stomach, and bladder and to a lower extent in brain, kidney, pancreas, and bone marrow. Expressed in mononuclear leukocytes and megakaryocytes.