

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TERF2IP Polyclonal Antibody

多克隆抗体

产品货号: PA05712

第三版

描述:The gene encodes a protein that is part of a complex involved in telomere length regulation. Pseudogenes are present on chromosomes 5 and 22.

抗原:Synthesized peptide derived from the Internal region of human RAP1.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous. Highly expressed.