

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat RASSF4 Polyclonal Antibody

多克隆抗体

产品货号: PA05714

第三版

描述:The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described. RASSF4 (Ras Association Domain Family Member 4) is a Protein Coding gene. Among its related pathways are Hippo signaling pathway -multiple species. An important paralog of this gene is RASSF2.

抗原:Synthesized peptide derived from the Internal region of human RASSF4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: