

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RASSF6 Polyclonal Antibody

多克隆抗体

产品货号: PA05715

第三版

描述: This gene encodes a member of the Ras-association domain family (RASSF). Members of this family form the core of a highly conserved tumor suppressor network, the Salvador-Warts-Hippo (SWH) pathway. The protein encoded by this gene is a Ras effector protein that induces apoptosis. A genomic region containing this gene has been linked to susceptibility to viral bronchiolitis. Alternative splicing results in multiple transcript variants and protein isoforms.

抗原: Synthesized peptide derived from the Internal region of human RASSF6

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highest expression in thymus, kidney and placenta. Also detected in colon, small intestine and lung. Tends to be down-regulated in 30-60% of tumors derived from breast, colon, kidney liver, rectum, pancreas, stomach and the thyroid gland compared to the normal counterpart.