

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat RXFP2 Polyclonal Antibody

多克隆抗体

产品货号: PA05730

第三版

描述:This gene encodes a member of the GPCR (G protein-coupled, 7-transmembrane receptor) family. Mutations in this gene are associated with cryptorchidism. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. RXFP2 (Relaxin/Insulin Like Family Peptide Receptor 2) is a Protein Coding gene. Diseases associated with RXFP2 include Cryptorchidism and Anorchia. Among its related pathways are Peptide ligand-binding receptors and CREB Pathway. GO annotations related to this gene include G-protein coupled receptor activity and protein-hormone receptor activity. An important paralog of this gene is RXFP1.

抗原:Synthesized peptide derived from the Internal region of human Relaxin Receptor 2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed mainly in the brain, kidney, muscle, testis, thyroid, uterus, peripheral blood cells and bone marrow.