

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CHML Polyclonal Antibody

多克隆抗体

产品货号: PA05731

第三版

描述:The product of the CHML gene supports geranylgeranylation of most Rab proteins and may substitute for REP-1 in tissues other than retina. CHML is localized close to the gene for Usher syndrome type II. CHML (CHM Like, Rab Escort Protein 2) is a Protein Coding gene. Diseases associated with CHML include Choroideremia. Among its related pathways are Vesicle-mediated transport and RAB geranylgeranylation. GO annotations related to this gene include GTPase activator activity and Rab geranylgeranyltransferase activity. An important paralog of this gene is CHM.

抗原:Synthesized peptide derived from the Internal region of human REP-2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: